
NON INVASIVE INFANT SLEEP APNEA DETECTION USER MANUAL

CO321 CO324 CO325 UNIFIED PROJECT

TEAM

E/14/158: GIHAN JAYATILAKA
E/14/339: SUREN SRITHARAN
E/14/379: HARSHANA WELIGAMPOLA
E/14/237: PANKAYARAJ PATHMANATHAN

*Deaprtment of Computer Engineering
Unnviersity of Peradeniya*

2018

Contents

1	Introduction	2
1.1	Sleep apnea	2
1.2	The system	2
1.3	Process overview	2
2	Features	2
3	Quick Start up	3
4	Device Installation	3
4.1	Power requirement	4
4.2	Network configuration	4
4.3	Configuring Communication	4
4.4	Configuring Camera	4
4.5	Device Reset	5
5	Starting the Node functionality	5
5.1	Starting services	5
6	Using the web interface	5
6.1	Registering device	6
6.2	Delete device	6
6.3	View Breathing pattern	6
7	Screenshots	7

List of Figures

1	Signup page	7
2	Login page	7
3	Breathing pattern reports	8
4	Notifications	8
5	Add device	9
6	List of devices	9
7	Device administration	10

1 Introduction

1.1 Sleep apnea

Sleep Apnea is a serious disorder caused by the interruption of breathing during sleep. This can cause the people to stop breathing for several times (even hundreds) if not treated properly. It can affect people of any age. But when the babies are affected with the condition they tend to not get up and keep on sleeping which may risk their lives.

1.2 The system

This is a non invasive solution based on video processing. The infant is observed by a video camera which is connected to a single board computer (Raspberry pi) which analyzes the video feed to diagnose breathing anomalies. The camera is turned to a proper orientation for the observation using a robotic arm.

1.3 Process overview

The system uses a high definition camera to monitor the baby and a powerful single board computer to process the video to identify the breathing pattern and asses it for anomalies in order to detect **sleep apnea** using intelligent algorithms. The unit communicates with the web server over the internet which provides an user friendly interface for the doctors and parents to see the reports.

2 Features

- Measurement and storage of breathing pattern
- Real time alerts
- Our solution is 100% non intrusive.
- The camera node could operate independently with the inbuilt battery cells. There is no need even for internet connection for the primary functionality.
- Our camera can automatically detect the baby and orient itself to the correct position intelligently using the robotic arm.
- The *breathing detection algorithm* is automated than the previous work found. (The existing algorithms require human intervention to detect interesting regions etc;)
- The *breathing detection algorithm* is accuracte than the existing algorithms.

3 Quick Start up

This section explains quicker way to test the main functionality of the system.

1. Website
 - (a) Go to <http://isad.teambitecode.com/>
 - (b) Signup
 - (c) Login
 - (d) Register the device.
2. Device
 - (a) Connect to device
 - (b) Open config.ini file.
 - (c) Specify user name, device id and access token (Same as registered device).
 - (d) Press reset button

Now the status LED must become stable after some time. If you press the analyze button next to the device in the device list menu you must be able to observe the graph being plotted.

4 Device Installation

- Connect the device to the computer LAN
- Use a **SSH client** to connect to the node
- Enter the following commands

```
wget http://teambitecode.com/projects/sleep-apnea-detection/node.tar
tar -xvf node.tar
cd node
export $SDCARD=[The mount location of the SD card here]
export $WIFI=[Wifi network ssid name here]
export $WIFIUSER=[Wifi username here]
export $WIFIPASS=[Wifi password here]
sudo chmod 700 install.sh
sudo ./install.sh
```

- Now the device could be used anywhere in the building via wireless network.

Note: In case the Wi-Fi router restarts, please reset the node as well. Make sure the router is functional and able to connect when the node starts.

4.1 Power requirement

The device requires to be powered using a 5v 2A (minimum) power supply. The OEM power supply (5v 2.5A) provided by RaspberryPi would be ideal, but any other power supply that fulfills the condition would suffice.

4.2 Network configuration

In case the network configuration changes (ie. the wifi network ssid or password changes), connect to the node via **ssh** and enter the following commands.

```
export $WIFI=[New Wifi network ssid name here]
export $WIFIUSER=[New Wifi username here]
export $WIFIPASS=[New Wifi password here]
sudo ./savewificonfig.sh
sudo reboot -h now
```

4.3 Configuring Communication

The system communicates with the web server, therefore if the communication must be successful, first the device must be registered.

1. Open config.ini in file.
2. Set username, device.id and access.token.

If the device has been registered in the website and the configurations have been set up properly, then the status LED must be steady after some time after resetting the device.

4.4 Configuring Camera

The camera can be controlled manually or through the infant detection algorithm + camera rotation mechanism provided. To chose between these two options follow these steps:

1. Open config.ini in file.
2. Change auto_cam value under cam.conf
 - 1 - automatic camera rotation
 - 0 - manual control

If the configuration has been done properly then the status LED must blink before becoming steady after pressing the reset button if the status LED blinks for some time before becoming steady, then the automatic camera rotation works fine. The status LED blinks during the rotation phase and after completion LED becomes steady.

Note: Please make sure that the camera is pointing towards the infant when the system starts up. If the infant is not in the line of sight of the camera, automatic rotation of the camera will not work as intended!

Note: The infant should now wear any cloths in the stomach area. Because this device uses the movement of the stomach to identify the breathing patterns. Therefore, make sure not to wear any cloths in abdomen (t-shirts) or put any cloths on top of the infant (towel or small pieces of cloth that are used to wrap around the infant)

Warning: Do NOT keep the device close to the infant. When the camera is too close, the video taken might not be suitable/out of focus to process. Also, if the device is too close, the infant might unknowingly tamper the device and the infant might get hurt.

4.5 Device Reset

The device reset button could be used in case the device malfunctions or if the configurations have been changed

1. Press the reset button.
2. Wait for the LED to turn on (blink/ steady).

5 Starting the Node functionality

5.1 Starting services

When the node boots to the Operating System root user, it automatically starts the services on system start-up. If this does not happens you can manually start the services by running the `services.py` python file on the default directory.

You can start the nodes functionality when the service is running. Press the button once to start the node. At the beginning the camera with rotate to focus on the infant. Afterwards, the algorithm will run while collecting breathing information.

Note: Do not change the camera position manually when the algorithm is running. If you want to change the camera position. Reset the device first.

6 Using the web interface

The device communicates with the server to provide live alerts and plot graphs. In order to facilitate this the user must sign up and the device must be registered.

1. Go to www.isad.teambitecode.com/

2. Signup by providing user name, email and password. (If you haven't created an account already)
3. Login using user name and password.

6.1 Registering device

Device registration requires a unique device ID and an access token which is used to authenticate the device. (Note that the access token cannot be changed (for security reasons), therefore if forgotten the device must be removed and registered back again).

1. Go to Add devices menu under Devices.
2. Enter Device ID (unique) and access token.
3. Click Register.

To ensure that the device has been registered go to the Device list menu under devices.

6.2 Delete device

1. Go to Device list menu under Devices.
2. Click delete button next to the device.
3. Click Ok to confirm.


To ensure that the device has been deleted go to Device list menu under devices.

6.3 View Breathing pattern

1. Go to Device list menu under Devices.
2. Click analyze button next to the device.


Now you should be able to view the graph being plot.

7 Screenshots



A screenshot of a web browser showing the 'create an account' page. The browser's address bar displays '127.0.0.1:8000/client/signup/'. The page has a teal background. A white form is centered, titled 'create an account'. It contains four input fields: 'What's your username?', 'Choose a password', 'Confirm password', and 'Email address'. Below these is a teal 'Sign me up!' button. At the bottom of the form, it says 'already have an account? [login](#)'.

Figure 1: Signup page



A screenshot of a web browser showing the 'account login' page. The browser's address bar displays '127.0.0.1:8000/client/login/'. The page has a teal background. A white form is centered, titled 'account login'. It contains two input fields: 'enter your username' and 'enter your password'. Below these is a teal 'Sign me in!' button. At the bottom of the form, it says 'create an account - [forgot password](#)'.

Figure 2: Login page




Figure 3: Breathing pattern reports




Figure 4: Notifications

Infant sleep apnea detection

Search...

Dashboard

Devices

Device List

Add Device

Patients

Patient chart

Devices

Add New Device

Device ID

Device ID must be unique for each device.

Access Token

Password

Register

Clear

Figure 5: Add device

Infant sleep apnea detection

Search...

Dashboard

Devices

Device List

Add Device

Patients

Patient chart

Devices

Device List

Show entries

Search:

Device No.	Device ID	Active	Last Update	Actions
1	dev1		-	Analyze Delete
2	dev2		-	Analyze Delete
3	dev3		-	Analyze Delete
4	dev4		-	Analyze Delete

Showing 1 to 4 of 4 entries

Previous1Next

Figure 6: List of devices

Django administration

Home : Authentication and Authorization : Users : Baby01

Change user

Username:

Baby01

Required. 150 characters or fewer. Letters, digits and @/./+/-/_ only.

Password:

algorithm: pbkdf2_sha256 iterations: 100000 salt: lKyes2***** hash: O+PZfD*****

Raw passwords are not stored, so there is no way to see this user's password, but you can change the password using this form.

Personal info

First name:

Baby

Last name:

One

Email address:

Baby1@gmail.com

Permissions

☒ Active

Designates whether this user should be treated as active. Unselect this instead of deleting accounts.

☐ Staff status

Designates whether the user can log into this admin site.

☐ Superuser status

Figure 7: Device administration

END